

PROMOTIONAL OFFERS*

Valid from 1st of September till 15th of December 2025

Only for participating dealers Benelux.
All prices are in €, exclusive VAT.



*Benelux area

Menu

OPUS	4
Silky	7
C20	10
Newtron	13
Piezotome Cube	17
Implant center	18
Miniled	20
Riskontrol	24
Mouth mirrors	26
Expasyl	27
X-mind Dc	28
X-mind Unity	29
Sopix 2	30
X-mind Optima 3D	32
X-mind Prime 2	36
PsPix	39
Contact	43

opus

Original air Polisher and Ultrasonic Scaler

€ 7.968,00



**Made in France.
Trusted worldwide.**

**OPUS® is fully engineered
in France,**

with our Made in France quality seal and our commitment to handmade precision.

Our French craftsmanship allows us to proudly master:

- Full design and development manufacturing
- Supply chain processes
- Customer care service

The choices empowering you to work the way you want

OPUS® device offers user-friendly experience with:

- Instant treatments to contribute to **efficient and safe operations** with its responsive color-guidance interface
- Simplified mode adjustments
- Power level regulation for various clinical indications



ColorDrive Control

Simplifies power selection for ultrasonic and polishing mode.

QuickSense

Recognize automatically handpiece and power selection.



Intuitive touch screen

Displays only available power per mode. Enhanced handling even with wet gloves.



Easy to keep clean

Smooth high-quality surface made of hygienic glass. Integrated tubes cleaning system.



Smart foot control

4 functions wired pedal with integrated handle.

Keep your practice running smoothly with single-piece pressurizer **EasyLock**:

- Eliminate frustrating clogs
- Ensure consistent powder flow
- No diver required
- Guarantee safety and instant polishing




Quick powder insertion and removal
Grab & Go mechanism for instant polishing and switch between powders.

You won't waste time or product with continuous pressure put into the bottle ensuring powders' residues to stay inside the tube and allow the usage of all the powder until the very last grain.

Our **single-use recyclable SILKY® powder bottles** eliminate the need for refills and problematic powder dust for faster powder handling with less downtime.

Our **single-use recyclable SILKY® powder bottles** eliminate the need for refills and problematic powder dust for faster powder handling with less downtime.



"The air polishing is a nice experience for a patient."

DH Carlotta Becciani (Italy)

Experience a new symphony with continuous flow...

OPUS® device is combining our legacy ultrasonics technology, reinvented air polisher with **clogging control by design** and new powders' offering to **provide versatility at your fingertips**.

Our new SILKY® powders embody **minimally invasive, gentle yet effective treatments**. This new range of powders is designed to provide an effective polishing solution for supra- and subgingival procedures.

SILKY® Supra



Sodium bicarbonate-based



Grain size: 40µm




Natural mint and lemon flavors for patients' comfort



Effectively remove biofilm and heavy stains

SILKY® Perio



Glycine-based powder



Grain size: 25µm



Gentle yet effective cleaning



Ensures gentleness on tissues while its low abrasiveness enhance safety

SILKY® pen: our lightest and slimmest handpiece ever made

As a dental professional, you strive for efficient, effective, and comfortable patient care. SILKY® pen, our new air-polisher designed to **effectively treat both sub- and supragingival indications** combined with its 4 dedicated SILKY® nozzles offering a unique combination of **versatility, precision, and ease of use.**

Featuring an **ultra-light and slim and clogging control by design**, the SILKY® pen handpiece provides a **secure and comfortable grip for enhanced control.**

“Having a device for air polishing can be of big added value in the practice.”

Dr. Joseph Sabbagh (Belgium)

SILKY® nozzles:

- Precise and continuous spray for supra- and subgingival biofilm removal
- Designed to control clogging and reduce interruptions and maintenance
- High-quality chrome cobalt made for durability and corrosion resistance

SUPRA Nozzle

120° wide angled spray for fast and efficient supra-gingival prophylaxis.

PERIO MAINTENANCE Nozzle

Refined 0.6mm powder outlet for gentle and precise periodontal maintenance.

PERIO EASY Nozzle

Innovative shoe-horn design for easy insertion and targeted cleaning of moderate pockets.

PERIO Nozzle

Automatic pressure reduction (<1 bar) for safe and effective deep pocket decontamination.

“First impressions were very positive... It performs outstandingly. Having the air polishing unit gives practitioners flexibility”

Dr. Rajan Nansi (UK)

ColorDrive system

The RFID tracking mechanism enables **automated powder recognition** with colored spotlight and pressure selection for SUPRA and PERIO modes allowing **effortless handling with complete peace of mind.**



C20

Easy to use. Easy to master.

More than just a tool, The C20 enhances the diagnostic process by creating an engaging and interactive experience for both the dental team and their patients.

By streamlining communication and boosting chairside efficiency, the C20 helps transform practices into high-performing, patient focused environments.



Join the game
OF MODERN
DENTISTRY



Level up
YOUR PRACTICE



Press start
AND GO

€ 1.199,00

MADE IN
EUROPE



Press start, And Go.

The C20 offers a **plug-and-play** experience, making it **straightforward** to install in any dental office. Its **intuitive design** ensures that users can start capturing and managing images immediately, enhancing workflow without any disruption.



Plug
& Play

Ready to use with
pre-installed driver.



Universal
connectors

USB 2 & universal
socket.



Customize
settings

A wide range of settings
to suit your needs.




Windows®
& Mac®
compatible



Exclusive and patented Color Coding System™

Power selection is accurately determined with the Color Coding System™, which intuitively associates each tip color ring with one of the 4 power ranges



newtron p5XS

1

Ref. F00246



1S

Ref. F00245



H3

Ref. F00369



P5XS b-LED incl. with the following tips:
Ref. F62100

~~€2788,25~~**€2099**

N°4
Ø22


Ref. MP3622PH

NOW INCLUDED


PURE REFLECT®
with a value
of €54,89

Ref. MP3622PH

Simple and intuitive settings with the Color Coding System™



Newtron® handpiece



ref. F12281 - Autoclavable

~~€635,50~~

€339

Newtron® LED handpiece



ref. F12609

~~€879,63~~

€495

Newtron® Slim B. LED handpiece




- with blue LED ring: ref. F12900
- with white LED ring: ref. F12905

~~€879,63~~

€549

hygiene

Hygiene Kit




ref. F00935

~~€529,14~~

€465

periodontics

Perio Kit




ref. F00936

~~€670,97~~

€589

implantprotect pure titanium

Implant Protect Kit



ref. F02120

~~€735,54~~

€669

endosuccess retreatment

Endo Success Retreatment Kit





ref. F00737

~~€702,36~~

€620

- Excellent view, full LED ring
- Avoid necrosis up to 120ml irrigation with peristaltic pump
- Improve durability of the handpiece, external irrigation
- Maximum hygiene, sterile single use or autoclavable lines
- Work don't play, use of big bottles (1l) possible
- Large choice of clinical application, 38 tips
- Tips II are compatible with IC2, Piezotome Cube, Piezotome2, Piezotome solo



piezotome cube

**Dynamic power adaptation for
safe bone surgeries**

~~€7.609,01~~
€5.650

Ref. F50100

Delivered with a multifunction footswitch, a Piezotome® CUBE LED handpiece, 1 handpiece holder, an Essential II kit (6 tips BS1S, BS4, SL1, SL2, SL3 et LC2), a Piezotome® dynamometric wrench, 1 irrigation pole, 2 autoclavable irrigation lines, 5 clips and 15 perforators.



Intraoral surgery ultrasonic control console.
Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

implant center2

3 modes:

- Piezotome® (bone surgery).
- Newtron® (perio, endo,...).
- I-Surge™ (implant motor).

Ref. F27200 - Delivered with a LED motor, a LED motor cord, a progressive footswitch, 2 brackets, 2 handpiece holders, 2 sterilizable irrigation lines and 30 perforators.

The alliance of technologies
for powerful and secure
bone surgeries



~~€14.538,36~~
€9.599

Ultrasonic control console for conventional dentistry and intraoral dentistry and motor control console for dental implantology
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

Accessoires

Promotion valable only for ImplantCenter2

Newtron® LED handpiece



Ref. F12609
For Piezotome® 2 and ImplantCenter™ 2

Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

Newtron® LED cord



Ref. F12830
For Piezotome® 2 and ImplantCenter™ 2



~~€2.022,97~~
€1.489

Universal irrigation lines Sterile x10



~~€184,07~~
€145

as from 4 pieces

Ref. F59905 - Single-use with integrated cassette and flow control system.

Class Is medical device (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France

PIEZOTOME® LED handpiece + cord



~~€1.635,16~~
€1.299

Ref. F12815

Class IIa medical devices (GMED) - CE 0459
Manufacturer: SATELEC® a company of ACTEON® group - France



Bone Surgery™ II Kit*

Ref. F87509 - Tips BS1S II, BS2L II, BS2R II, BS4 II, BS5 II, BS6 II.

~~€1.432,35~~ **€1.090**



Sinus Lift II Kit*

Ref. F87519 - Tips SL1 II, SL2 II, SL3 II, SL4 II, SL5 II.

~~€1.077,92~~ **€990**



Universal II Kit*

Ref. F87539 - Tips BS6, LC1, LC2, CE1, CE3 et PZ1.

~~€1.393,47~~ **€1.090**



Essential II Kit*

Ref. F87528 - Tips BS1S II, BS4 II, SL1 II, SL2 II, SL3 II, LC2 II.

~~€1.432,35~~ **€1.090**



Intralift™ II Kit*

Ref. F87536 - Tips TKW1 II, TKW2 II, TKW3 II, TKW4 II, TKW5 II (Ø 3 mm).

~~€1.367,71~~ **€1.090**



Extraction™ II Kit*

Ref. F87546 - Tips LC1 II, LC2 II x2, LC2L II, LC2R II, Ninja® II.

~~€1.087,60~~ **€990**



Crest Splitting II Kit*

Ref. F87567 - Tips CS1 II, CS2 II, CS3 II, CS4 II, CS5 II, CS6 II.

~~€927,18~~ **€819**



Crown Extension II Kit*

Ref. F87554 - Tips BS6 II, CE1 II, CE2 II, CE3 II.

~~€737,70~~ **€619**

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.
Class IIa medical devices (GMED) - CE 0459 - Manufacturer: SATELEC® a company of ACTEON® group - France

* Kits are delivered with an autoclavable metal tip holder and an autoclavable universal wrench.

miniled ortho 2

€1.450,04

€1.090

Ref. F05220

• Customizable signal

An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec.

• Power

3000mW/cm²
(Ø5,5mm light guide).

• Long lasting battery

MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day.

• Ortho friendly

Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch.

• Efficiency & safety

Cure any type of material used orthodontics efficiently and quickly.

No noisy and unhygienic fan needed to cool the lamp.

Specialized Orthodontic Treatments

Cure in a simple gesture, in few seconds from few brackets to a whole quadrant



Simply Powerful

Outperforms in
power, efficiency and speed

€1.400,53

€1.090

Ref. F05214



miniled supercharged

• Visibility

Luminous screen with customizable modes and times, count curing time and load indicator.



• Power

2000mW/cm²
(Ø7,5mm light guide).

3000mW/cm²
(Ø5,5mm light guide - optional).


• Optimal result

3 modes to fit all clinical needs.




FAST MODE

The choice of rapidity
Full power during 3, 4, 5 or 10 sec, according to the composite thickness.



PULSE MODE

Limits temperature rise
5 or 10 flashes of 1 sec for a polymerization in security.



RAMPING MODE

Avoids composites retraction
Progressive power increase during 9 or 15 sec to ensure perfect sealing of restoration.

miniled standard

• Fast

Only 6 to 12 sec to polymerize
2mm composite.

• Power

1250mW/cm²
(Ø7,5mm light guide).


2000mW/cm²
(Ø5,5mm light guide - optional).

• Safety

Base-charger includes a radiometer for light power
intensity check-up.

• Optimal result


3 modes to fit all clinical needs.



FAST MODE

Fast and easy


Full power during 10 sec
(bip at 5 sec) on regular
depth.



PULSE MODE

Ally rapidity & safety

10 flashes of 1 sec
(bip after 5 flashes)
with times off to
preserve dental pulp.



RAMPING MODE

No stress cracking

Progressive power
increase during 20 sec
(bip every 5 sec) for
a complete and safe
curing.

~~€1.090,24~~

€890

Ref. F02507

Highly Adaptable

The alliance of power & safety



servotome

Delivered with a single footswitch,
a box with 10 electrodes,
a conductive bracelet.

Ref. F57210

~~€2.025,86~~

€1.699



Riskontrol® contributes to clinical protection



RISKONTROL® is also a GREEN solution

- Disposable tip
- Economic: simplify process to save time
- Recyclable plastic: environmentally friendly
- Cost-effective
- Flexible for a better use: respect of patient's anatomy
- Safety lock system: no risk of ejecting the tips




Sterilizable tips are not sufficient anymore

RISKONTROL® can be used for all procedures in any workflow




Easy to fit

1




It's installed.

2



It's fitted.

3



It's locked.

Discover our RISKONTROL® range

Adapt your tips to the taste
of your patients



by **art riskontrol** ~~€61,53~~
€32,50

- Ref. 201 700 - Mint
- Ref. 201 701 - Mandarin
- Ref. 201 702 - Liquorice
- Ref. 201 703 - Anise
- Ref. 201 704 - Black currant

Non-sterile disposable tips (x 250)

Protect entirely your air/water syringe
thanks to the plastic sheath



total protec
by **riskontrol**

~~€32,93~~
€29,50

- Ref. 201 750
- Ref. 201 751

Non-sterile disposable tips with plastic sheath (x 100)

**When you buy
10 boxes 1 adaptor FOC**

Add color

to your patients' care experience



~~€56,65~~
€32,50

by **classic riskontrol**

- Ref. 201 664 - White
- Ref. 201 665 - Blue
- Ref. 201 666 - Yellow
- Ref. 201 667 - Green
- Ref. 201 668 - Pink

Non-sterile disposable tips (x 250)

Choose a totally sterile
workflow



~~€38,13~~
€30,50

sterile
by **riskontrol**

- Ref. 201 770

Sterile disposable tips (x 96)

Classic by RISKONTROL®, Art by RISKONTROL®, Total Protec by RISKONTROL®: Class I medical device. Sterile by RISKONTROL®: Class Is medical device. CE0459 (GMED). Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of these devices is in the Product Feature Summary available on the laboratory's website. Manufacturer: PRODUITS DENTAIRE PIERRE ROLLAND - ACTEON® Group - FRANCE.

Mouth mirrors

PURE REFLECT®

- Very bright image
- Resistant

For all purposes

RH

- A greyish reflection



TI

- A blueish reflection

FRONT SURFACE

Cross-sectional view






Image reflected directly on the surface of the glass, distortion-free.
NO DOUBLE IMAGE

BACK SURFACE

Cross-sectional view of an aluminium mirror



Reflecting surface under the glass. **DOUBLE IMAGE**



PURE REFLECT®

N°	Ref.	N.
4	MP3622PH	12
5	MP3624PH	12

~~€54,89~~
€39

Mirrors: Class I professional dental medical devices - CE
Manufacturer: PRODONT HOLLIGER ACTEON® Group - France

A TRIPLE ACTION



Kaolin clay

1 ASTRINGENT EFFECT

helps to completely dry the sulcus with no blood or fluids.



2 OPEN THE SULCUS

effectively with no pain or damage to tissues thanks to the viscosity carefully calculated with our clinical expertise.



Aluminum chloride

3 PRESERVE EPITHELIAL ATTACHMENT*

for a reduced risk of gingival recession and bone resorption.

*Canula must not be inserted into the sulcus

expasyl
exact



~~€144,33~~
€99

Paste for temporary sulcular opening
20 capsules (0,3g) Expasyl™ Exact - Ref. 261010



expasyl

~~€250,09~~
€139

Paste for temporary sulcular opening
20 capsules (1g) Expasyl™ - Ref. 261030



~~€111,70~~
€79

Straight - Ref. 261005
Bended - Ref. 261015

Class I medical devices - CE - Manufacturer: PRODUITS DENTAIRE PIERRE ROLLAND - ACTEON® Group - France

x-mind dc



The X-MIND® generators can accomodate any operatory configuration.

3 arms lengths are available : 41 cm (16.2''), 82.5 cm (32.5'') and 110 cm (43.5'') and can be:

- Top wall mounted
- Bottom wall mounted

Dental unit, ceiling and mobile version are also available.

With its adaptable wall plate, it is easy to replace an older generator with the X-MIND® whitout having to do additional drilling on the wall.

~~€4.611,95~~
€3.150

Available with Remote Wireless Control



No more wiring issues

Say goodbye to the complicated and expensive wiring work traditionnally needed to connect the remote button in your practice.



Cost effective

Save on installation expenses with the easy-to-integrate wireless solution.



Enhanced convenience

Enjoy the freedom and flexibility of a wireless connection, making your practice more efficient and user-friendly.

Wireless Button Kit X-MIND DC

Ref. 5.15.00135

All X-Mind DC devices delivered from April 2021 onwards

€249



Wall version	Top mounting	Top mounting square cone	Bottom mounting	Bottom mounting square cone
0,40 m extensor	Ref. 5.92.23137	Ref. 5.92.23167	Ref. 5.92.23130	Ref. 5.92.23160
0,80 m extensor	Ref. 5.92.23138	Ref. 5.92.23168	Ref. 5.92.23131	Ref. 5.92.23161
1,10 m extensor	Ref. 5.92.23139	Ref. 5.92.23169	Ref. 5.92.23132	Ref. 5.92.23162

All references for the X-MIND® DC

Choose from a variety of cones to better fit to your clinical needs



Long cone (31 cm - 12'')
Recommended for the parallel technique



Short cone (20 cm - 8'')
Recommended for bisector technique



Rectangular Long cone (31 cm - 12'')
Reduces the X-rays applied to the patient's skin by 50%.

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy

x-mind unity

~~€5.767,06~~
€3.990

A sharp and contrasted image

- A contrasted and defined image
- Optimized protection for the patient
- The generator and sensor association

Fixing from above, arm 0.80 m, short cone.
Other configurations: see public price 2024.

Intraoral X-ray

- Very high quality
- Precise and flexible positioning
- Fast selection from traditional to digital mode
- Minimum exposure time: 0,02s
- Programmable timer
- Top or bottom mounting

Class IIb medical devices (IMQ) - CE0051 Manufacturer: de Götzen ACTEON® Group - Italy





sopix 2



A remarkable image quality, even for endodontics and the more complexes

SOPIX^{®2} provides a high resolution image by exploiting the latest technologies available on the market.



SOPIX^{®2} USB Sensor size 1
high definition with a sensor
holder, 10 sensor sheaths
Ref. S_802_0008

~~€6.920,49~~ **€3.990**

Also available in size 2
Ref. S_802_0009

~~€8.181,45~~ **€5.490**

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON[®] Group- France

Free yourself from constraints
with the embedded system and its
detachable sensor part.



SOPIX^{®2} Plug-in kit S1

Ref. SK60-1000

~~€6.524,54~~ **€5.990**


SOPIX^{®2} Plug-in kit S2

Ref. SK60-1001

~~€8.580~~ **€7.490**

Size 2
31x42mm

Size 1
25x39mm



sopix 2 plug-in

The latest technology with no
compromise on image quality with
more safety, freedom and comfort.

Resolution: >18lp/mm

Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON[®] Group- France

x-mind optima 3D

The Smartest, The Smallest, The Ultimate.

In today's rapidly evolving field of dentistry, the integration of CBCT and Artificial Intelligence is transforming patient care, enabling dental professionals to make more accurate diagnoses and perform safer procedures.

Embracing this cutting-edge integration is not only elevating the standard of care but also optimizing efficiency and patient communication in modern dental practices.

X-MIND® optima 3D represents the perfect combination of advanced hardware technology, enclosed in a lightweight and compact design and next-generation software, RealGuide™ 5.4, bringing the benefits of Artificial Intelligence into every step of patient care, from a precise diagnosis to more effective treatments and patient communication.

Meet the X-MIND® optima 3D Family.

—€81.999—
€44.499

* Option:
The additional cost for
the 17x12 version is €3,000

X-MIND® optima 2D/3D wall 100-240V

2D/3D PAN device 5x5 to 12x10 wall installation
Ref. W1203005



—€93.499—
€54.999

* Option:
The additional cost for
the 17x12 version is €3,000

X-MIND® optima 2D/3D ceph 100-240V

2D/3D PAN CEPH device 5x5 to 12x10 floor installation
Ref. W1203010




—€83.499—
€45.999

* Option:
The additional cost for
the 17x12 version is €3,000

X-MIND® optima 2D/3D floor 100-240V

2D/3D PAN device 5x5 to 12x10 floor installation
Ref. W1203005.FL

New



Exceptional clarity.

At the heart of the X-MIND® optima 3D is **state-of-the-art detector technology**, which guarantees exceptional image clarity and outstanding anatomical coverage. Its **comprehensive range of Fields of View (FOV)**, from 5x5 to 17x12* cm (Ø x H), makes it suitable for a **wide range of dental procedures**, from general dentistry to orthodontics and maxillofacial surgery, while increasing tolerance for patient positioning and reducing the risk of marginal artifacts.




EXTRAORDINARY CLARITY

Get the most out of quality with the brand-new **ClearIQ** filter



OUTSTANDING VERSATILITY

From 5x5 to 17x12*, a tailored solution for every clinical need



TOP-LEVEL PERFORMANCE

Enhanced image depth with the new kV output range 60-90




MAXIMUM PATIENT COMFORT

Quick and easy acquisition workflow




AUTOMATIC POSITIONERS RECOGNITION

Maximum accuracy



ZERO FOOTPRINT


Maximum accessibility



* Option - The additional cost for the 17x12 version is €3,000


5x5 cm FOV

- Improves the **detection of root canals**, particularly in case of abnormal root canal shapes



8.5x5 cm FOV

- Useful for **implant surgery** and in **general dentistry** as a low dose exam as limited to a single jaw




8.5x9 cm FOV

- Useful for **complex implant surgery** as well as left or right TMJ diagnosis



12x10 cm FOV

- Suitable for most **oral surgery cases** including multiple implant surgery



17x12* cm FOV

New

- Represents the perfect solution for **complex orthognathic cases**, **zygomatic implant planning** and **orthodontic surgery**



* Option - The additional cost for the 17x12 version is €3,000

Meet the 2D Prime Family

X-MIND® prime is available in three configurations: 2D wall mounted, 2D Ceph and 2D floor to suit all installation and clinical demands.



€13.999,00

X-MIND® prime WALL 110-240V

Order code: W1200001

The ACTEON®s smallest extraoral 2D imaging solution ever for the maximum **accessibility in any context.***



with one sensor: €24.999,00

with two sensors: €29.999,00

X-MIND® prime CEPH 110-240V

Order code: W1200009

Cephalometric add-on for a full 2D system still encompassing the same philosophy of prime, **ideal for pediatric applications.**



€14.999,00

X-MIND® prime FLOOR 100-240V

Order code: W1200001.FL

*Belgian market only (W1200001/W120001.PB)

Rely on a comprehensive set of 2D exams


X-MIND® prime provides a comprehensive range of 2D exams designed to meet the diverse needs of general dentistry, including dental panoramic imaging, temporomandibular joint assessments and maxillary sinus evaluations.



Dental panoramic



Maxillary sinuses



Temporomandibular joint



Child panoramic



Half-panoramic
Left & right



Bitewing
Left & right



Frontal

PSPiX²

is AUTONOMOUS


The fully automated process free yourself from constraints

Just drop the plate inside the slit, the PSPiX² will automatically detect the size, scan the plate, optimize the image, and send it instantly to your imaging software.

In the meantime, it will automatically erase and eject the plate for an immediate reuse or switch to standby mode.

**Scan
simply &
quickly**

Fully automated process



—€8.358,04—

€6.490

Ref. S_702_0001

1

Auto-access:

The door only opens when the imaging plate is inserted.

2

Auto-detect:

PSPiX² automatically detects the size of the inserted imaging plate.

3

Auto-scan:

PSPiX² scans and optimizes the image before displaying it into the right workstation.

4

Auto-eject:

The imaging plate is automatically erased and ejected for an immediate reuse.

5

Eco-mode:

PSPiX² automatically switches to economic mode when not in use for an extended period of time.

Protective bag for ERLM*

For **PSPIX®2** and **PSPIX®**
(previous generation)

~~€84,14~~
€57,40



	PSPIX® (previous generation)	
	Box of 250	
Size 0		700338L
Size 1		700339L
Size 2		700340L
Size 3		700341L

Protections for ERLM*

Protective cover for
PSPIX® and **PSPIX®2**
(previous generation)

~~€47,12~~
€26,25



Compatible with	PSPIX2	PSPIX
	✓	✓
PSPIX®	✓	✓
SOREDEX® Digora	✓	✓
INSTRUMENTARIUM® Express	✓	✓
DÜRR DENTAL® Vistascan®	✓	✓
KAVO®	✓	✓



	PSPIX® (previous generation)	
	Box of 250	
Size 0		700342L
Size 1		700343L
Size 2		700344L
Size 3		700345L

ERLM plates*

For **PSPIX®2**

~~€131,80~~
€95



	PSPIX®	
	Box of 2	
Size 0 (22 x 35mm)		990215
Size 1 (24 x 40mm)		990216
Size 2 (31 x 41mm)		990217
Size 3 (27 x 54mm)		990218

Protective bag for ERLM*

For **PSPIX®2**

~~€98,88~~
€70



	PSPIX®	
	Box of 300	
Size 0		700546
Size 1		700547
Size 2		700548
Size 3		700549

- Protection against scratches, light and cross contamination
- Ergonomic and easy to use
- Economic solution

* Luminescent Radio Memory Display for digital phosphor plate systems


PSPIX® imaging plate - Class IIa Medical devices (GMED) - CE0459
Manufacturer: Sopro ACTEON® Group - France
Bag & cover: Class I medical devices - CE
Manufacturer: Sopro ACTEON® Group - France



The thin and flexible imaging plates,
available in 4 sizes, will allow you to
easily scan all morphologies...



... and all clinical cases



CONTACTS

Area Sales Manager Equipment BeLux

Raphaël VERLINDE
+32 (0) 494 116 499
raphael.verlinde@acteongroup.com

Area Sales Manager Equipment The Netherlands

Roger VAN EIJDEN
+31 (0) 624 835 367
roger.vaneijdsen@acteongroup.com

Area Sales Manager Imaging BeNeLux

Andy LUBERTI
+31 (0) 613 948 984
andy.luberti@acteongroup.com

Imaging Specialist BeNeLux

Pierre DEKUSSCHE
+32 494 111 253
pierre.dekussche@acteongroup.com

Administration Commerciale

Sabine PARIS
Belgium
+33 (0) 556 349 332
Fax: +33 (0) 566 349 331
sabine.paris@acteongroup.com

Nirina INISAN
Luxembourg / The Netherlands
Tel.: +33 (0) 556 349 250
Fax: +33 (0) 566 349 331
nirina.inisan@acteongroup.com



ACTEON® CUSTOMER SERVICE

ACTEON's* customer service
can be contacted at:

Imaging
supportbenelux@acteongroup.com

Equipment / Pharma
cs@acteongroup.com
supportbenelux@acteongroup.com

ACTEON®
17 avenue Gustave Eiffel
33700 MERIGNAC - FRANCE
TEL: +33 (0) 556 340 607
FAX: +33 (0) 556 349 292
info@acteongroup.com
www.acteongroup.com

Non contractual document
Copyright© 2025 ACTEON®. All rights reserved. No information or
part of this document may be reproduced or disclosed in any format
without the prior authorisation of ACTEON®.
For professional dental use only.
Read the instructions in the notice carefully before use.
All information that is essential to ensure correct use of these
devices is included in the user manual or the summary of the
product characteristics available on www.acteongroup.com.
Creation date: January 2025

x-mind optima3D

EASY-JOIN QR Code

SCAN ME!



 **acteon club**